

# Dental LT Clear

Class IIa Long-Term Biocompatible Resin for Form 2 [\\$399 / L](#)

Formlabs Dental LT Clear Resin produces strong, accurate, biocompatible parts ideal for dental splints, retainers, and similar applications, and is specifically designed to work with your Form 2 3D printer.



FLDLCL01

**formlabs** 

**Prepared** 10 . 04 . 2017  
**Rev** 01 10 . 04 . 2017

To the best of our knowledge the information contained herein is accurate. However, Formlabs, Inc. makes no warranty, expressed or implied, regarding the accuracy of these results to be obtained from the use thereof.

# Material Properties Data

	METRIC	METHOD
	Postcured	
<b>Flexural Properties</b>		
Ultimate Flexural Strength	≥ 50 MPa (no break)	ISO 20795-2:2013
Flexural Modulus	≥ 1300 Mpa	ISO 20795-2:2013
<b>Hardness Properties</b>		
Hardness Shore D	80 - 90D	ISO 868:2003
<b>Impact Properties</b>		
Maximum stress intensity factory	≥ 11 MPa·m <sup>1/2</sup>	ISO 179:2010
Total fracture work	≥ 250 J/m <sup>2</sup>	ISO 20795-2:2013

Dental LT Clear is tested at NAMSA, Chasse sur Rhône in France, and is certified biocompatible per EN-ISO 10993-1:2009/AC:2010. Further details are available upon request.

The product is in compliance with ISO Standards:

- EEN ISO 1641:2009
- EN-ISO 10993-1:2009/AC:2010
- EN-ISO 10993-3:2009
- EN-ISO 10993-5:2009
- EN 908:2008